

VENOSAN®

LEGLINE 20

Support Stockings 140 den

ELEGANT LUXURY



VENOSAN®

LEGLINE 20

20 mm Hg
140 den

Support Stockings 140 den

HEEL AND TOE REINFORCEMENT



Stützstrümpfe | Support stockings | Bas de support | Calze riposanti

AD
Wadenstrümpfe
Below knee
Bas au mollet
Gambaletti

AGH
Schenkelstrümpfe
Thigh length
Bas à cuisse
Calza alla coscia

AT*
Strumpfhose
Pantyhose/Tights
Collants
Collant

ATU Materna*
Schwangerschaftsstrumpfhose
Maternity pantyhose
Collants de grossesse
Collant gravidanzaa

* **Komfort-Hosenteil**
Comfort panty
Culotte confort
Corpino comfort

Stufenlos abnehmende Kompression | Graduated compression |
Compression dégressive | Compressione graduata

Umfang | Circumference | Circonférence | Circonferenza (cm)

AD	AGH, AT, ATU	
	60 % 12 mm Hg	cF
80 % 16 mm Hg	80 % 16 mm Hg	cD
		cC
100 % 20 mm Hg	100 % 20 mm Hg	cB

XS	S	M	L	XL
42-44	45-47	48-50	51-53	54-60
29-31	32-34	35-37	38-40	41-43
31-33	34-36	37-39	40-42	43-45
20-21	22-23	24-25	26-27	28-29

Effective Compression

Verhindert geschwollene Beine
Prevents swollen legs
Évite des jambes enflées
Impedisce le gambe gonfie

Everyday Elegance

Transparente Optik
Transparent look
Look transparent
Look trasparente

Optimal Comfort

Für höchsten Tragekomfort, rutscht nicht
For maximum comfort, does not slip
Pour un confort optimal, ne glisse pas
Non stringe, non scivola

Spitzenhaftband | Lace top band | Adhésive ornée | Bordo autoreggente pizzo




Aussen
Outside
Face externe
Lato esterno



Innen
Inside
Côte intérieur
Lato interno

Grösse | Height | Grandeur | Altezza



Artikel | Article | Article | Articolo

Type	*Colours	XS	S	M	L	XL
AD	Nude	VL 303 NU	VL 304 NU	VL 305 NU	VL 306 NU	-
	Sahara	VL 303 SA	VL 304 SA	VL 305 SA	VL 306 SA	-
	Black	VL 303 BL	VL 304 BL	VL 305 BL	VL 306 BL	-
AGH	Nude	VL 323 NU	VL 324 NU	VL 325 NU	VL 326 NU	-
	Sahara	VL 323 SA	VL 324 SA	VL 325 SA	VL 326 SA	-
	Black	VL 323 BL	VL 324 BL	VL 325 BL	VL 326 BL	-
AT	Nude	VL 333 NU	VL 334 NU	VL 335 NU	VL 336 NU	VL 337 NU
	Sahara	VL 333 SA	VL 334 SA	VL 335 SA	VL 336 SA	VL 337 SA
	Black	VL 333 BL	VL 334 BL	VL 335 BL	VL 336 BL	VL 337 BL
ATU	Sahara	VL 343 SA	VL 344 SA	VL 345 SA	VL 346 SA	VL347 SA
	Black	VL 343 BL	VL 344 BL	VL 345 BL	VL 346 BL	VL 347 BL

*Colours



Nude



Sahara



Black



Materialzusammensetzung | Composition | Composition | Composizione

78 % Polyamide | 22 % LYCRA® Elastane





VENOSAN®
Compression Stockings

www.venosan.com

VENOSAN® is a registered
trademark of SWISSLASTIC AG
ST. GALLEN, Switzerland

an  Company



SWISSLASTIC AG ST. GALLEN
Rorschacher Strasse 304
CH-9016 St. Gallen
Switzerland