


VENOSAN®

5000

Compression Stockings and Armsleeves

THE MAGIC OF SILVER



Swiss Made

VENOSAN®

5000


VENOSAN® 5001
Ccl.I
18–21 mm Hg
(2,4–2,8 kPa)

VENOSAN® 5002
Ccl.II
23–32 mm Hg
(3,1–4,3 kPa)

VENOSAN® 5003
Ccl.III
34–46 mm Hg
(4,5–6,1 kPa)

Compression Stockings and Armsleeves

The premium product for those with higher expectations.



Pure silver

Contains silver-coated X-Static® fibre.
No reduction of the properties of silver even after 150 washing cycles.

Heat conductive

Silver is well known for its excellent heat conductive properties. In combination with the microfibre TACTEL®, silver produces rapid dissipation of skin moisture.

Durable


Suitable for active men and women, thanks to the innovative yarn combination of silver and cotton.

Colours

Due to technical printing reasons the colours may vary.



Black



Sand



Marokko










Mexico



40 % TACTEL® Polyamide
30 % Elastane
20 % Cotton/Baumwolle
10 % X-Static® – The Silver Fiber™




TACTEL® fibre is a registered trademark of INVISTA and is used under license to NILIT® for nylon fibre products.

RAL Gütezeichen not valid for armsleeves!

VENOSAN® 5000 Compression Stockings	COMPRESSION CLASS	COLOUR	AD 	ADH 	AG 	AGH 	AGG 	AT 	ATU 	
	Ccl.I	mexico		●	●	○	●	●	●	○
		marokko		▲	▲		▲		▲	
		sand		▲			▲		▲	
		black		●	●	○	●	○	●	○
	Ccl.II	mexico		●	●	○	●	●	●	○
		marokko		●	●	○	●		●	○
		sand		●	○	○	●		●	○
		black		●	●	○	●	○	●	○
	Ccl.III	mexico		●	○	○	○	○	○	○
marokko			▲							
black			▲							

VENOSAN® 5000 Compression Armsleeves	COMPRESSION CLASS	COLOUR	CG 	CG-H 
	Ccl.II	mexico		○

- OTS (off the shelf) and MTM (made to measure)
- only MTM (made to measure)
- ▲ only OTS (off the shelf)



Inner side of VENOSAN® 5000 under a microscope

HOW SHOULD YOU TAKE MEASUREMENTS?


Compression Stockings

For correct measurements it is very important to take your time and proceed systematically. Whenever possible, the measurements should be taken when the leg is rested and not swollen. The measurements should always be taken with the patient standing. Well fitting medical compression stockings are the key to a successful treatment and satisfied patients.

MEASURE:



AD
Below knee



ADH
Below knee
with self
supporting
plain top

Below knee AD, ADH

1. the circumference of the ankle cB
2. the circumference of the calf cC
3. the circumference of the leg just below the knee cD
4. the length of the leg from the base of the heel to just below the knee lD



AGH
Thigh length
with self
supporting
lace top



AGH
Thigh length
with self
supporting
plain top

Thigh length with self supporting top AGH, Thigh length with belt left or right AGG

1. the circumference of the ankle cB
2. the circumference of the calf cC
3. the circumference of the leg just below the knee cD
4. the circumference of the thigh cG
5. the inside leg lG (for AGG up to the groin lK)



AGG
Thigh length
with belt



AT
Pantyhose


Thigh length for Pantyhose AT

1. the circumference of the ankle cB
2. the circumference of the calf cC
3. the circumference of the leg just below the knee cD
4. the circumference of the thigh cG
5. the length up to the groin lK



Measurement Position / Circumference Measurements

- cG** Thigh (5 cm below groin) ————— **cG**
- cD** Narrowest point approx. 2–3 cm below the knee ————— **cD**
- cC** Widest part of the calf ————— **cC**
- cB** Ankle (narrowest point above ankle bone) ————— **cB**



Length Measurements

- lK** Up to the groin
- lG** The inside leg
- lD** The length of the leg from the base of the heel to just below the knee

CIRCUMFERENCE MEASUREMENTS (CM)

	SMALL	MEDIUM	LARGE	X-LARGE
cB Ankle	20–23	23–26	26–29	29–31
cC Mid-calf	29–39	33–42	36–45	39–47
cD Below the knee	28–38	32–41	35–44	38–46
cG Thigh	48–63	54–67	58–72	63–74

LENGTH MEASUREMENTS (CM)

AD/ADH	lD	short	35–40
		long	37–42
AGH	lG	short	63–72
		long	68–78
AGG	lK	short	67–74
		long	74–85
AT	lK	short	67–74
		long	74–85

SIZE

It is important that we achieve the correct pressure gradient in accordance with the compression class prescribed by the doctor. The pressure exerted on the leg is determined by 2 factors:

- a) Physical law (LaPlace), which says that with constant tension the pressure is reduced with increased circumference and vice versa.
- b) With increased stretch the tension is increased and vice versa. One effect is therefore partially compensated through the other: Larger circumference = pressure loss due to circumference, but pressure gain due to increased tension of the material.

Thanks to more than 50 years experience we are familiar with these effects and can therefore offer a simple size range in accordance with our sizing chart, which guarantees medical efficacy.

For legs with a special shape or form we offer a made-to-measure service with extra fast delivery. Please ask for the special order form for made-to-measure stockings (REF 995220E).

LENGTH

VENOSAN® 5000 has a very high stretch in the length and adapts to almost every leg length without slipping.

Subsequently compare with the size chart for:


1. Which size does the ankle circumference cB fall into (small, medium, large or x-large)?
2. Which size do the other measurements fall into cC, cD and cG?

Custom made stockings may not be returned for credit. The manufacturer cannot be held responsible for incorrect measurements.

5001

Compression Stockings Ccl.1

VENOSAN® 5001
Ccl.1
18–21 mm Hg
(2,4–2,8 kPa)






ARTICLE	COLOUR	LENGTH	TOE	SELF SUPPORTING TOP	S	M	L	XL	
AD	mexico	long	open		SG 45001	SG 45002	SG 45003	SG 45004	
		short			SG 45021	SG 45022	SG 45023	SG 45024	
		long	closed		SG 45011	SG 45012	SG 45013	SG 45014	
		short			SG 45031	SG 45032	SG 45033	SG 45034	
	marokko	long	open		SM 45001	SM 45002	SM 45003	SM 45004	
		short			SM 45021	SM 45022	SM 45023	SM 45024	
		long	closed		SM 45011	SM 45012	SM 45013	SM 45014	
		short			SM 45031	SM 45032	SM 45033	SM 45034	
	sand	long	open		SD 45001	SD 45002	SD 45003	SD 45004	
		short			SD 45021	SD 45022	SD 45023	SD 45024	
		long	closed		SD 45011	SD 45012	SD 45013	SD 45014	
		short			SD 45031	SD 45032	SD 45033	SD 45034	
	black	long	open		SN 45001	SN 45002	SN 45003	SN 45004	
		short			SN 45021	SN 45022	SN 45023	SN 45024	
		long	closed		SN 45011	SN 45012	SN 45013	SN 45014	
		short			SN 45031	SN 45032	SN 45033	SN 45034	
	ADH	mexico	long	open	plain	SG 45041	SG 45042	SG 45043	SG 45044
			short			SG 45061	SG 45062	SG 45063	SG 45064
			long	closed	plain	SG 45051	SG 45052	SG 45053	SG 45054
			short			SG 45071	SG 45072	SG 45073	SG 45074
		marokko	long	open	plain	SM 45041	SM 45042	SM 45043	SM 45044
			short			SM 45061	SM 45062	SM 45063	SM 45064
			long	closed	plain	SM 45051	SM 45052	SM 45053	SM 45054
			short			SM 45071	SM 45072	SM 45073	SM 45074
black		long	open	plain	SN 45041	SN 45042	SN 45043	SN 45044	
		short			SN 45061	SN 45062	SN 45063	SN 45064	
		long	closed	plain	SN 45051	SN 45052	SN 45053	SN 45054	
		short			SN 45071	SN 45072	SN 45073	SN 45074	

ARTICLE	COLOUR	LENGTH	TOE	SELF SUPPORTING TOP	S	M	L	XL
AGH	mexico	long	open	plain	SG 45611	SG 45612	SG 45613	SG 45614
		short			SG 45641	SG 45642	SG 45643	SG 45644
		long	closed	lace	SG 45601	SG 45602	SG 45603	SG 45604
		short			SG 45631	SG 45632	SG 45633	SG 45634
	marokko	long	open	plain	SM 45611	SM 45612	SM 45613	SM 45614
		short			SM 45641	SM 45642	SM 45643	SM 45644
		long	closed	lace	SM 45601	SM 45602	SM 45603	SM 45604
		short			SM 45631	SM 45632	SM 45633	SM 45634
	sand	long	open	plain	SD 45611	SD 45612	SD 45613	SD 45614
		short			SD 45641	SD 45642	SD 45643	SD 45644
		long	closed	lace	SD 45601	SD 45602	SD 45603	SD 45604
		short			SD 45631	SD 45632	SD 45633	SD 45634
	black	long	open	plain	SN 45611	SN 45612	SN 45613	SN 45614
		short			SN 45641	SN 45642	SN 45643	SN 45644
		long	closed	lace	SN 45601	SN 45602	SN 45603	SN 45604
		short			SN 45631	SN 45632	SN 45633	SN 45634
	AGG left	mexico	long	open	SG 45151	SG 45152	SG 45153	SG 45154
	AGG right		long		SG 45171	SG 45172	SG 45173	SG 45174
	AT	mexico	long	closed	SG 45221	SG 45222	SG 45223	SG 45224
			short		SG 45231	SG 45232	SG 45233	SG 45234
		marokko	long	closed	SM 45221	SM 45222	SM 45223	SM 45224
			short		SM 45231	SM 45232	SM 45233	SM 45234
		sand	long	closed	SD 45221	SD 45222	SD 45223	SD 45224
			short		SD 45231	SD 45232	SD 45233	SD 45234
black		long	closed	SN 45221	SN 45222	SN 45223	SN 45224	
		short		SN 45231	SN 45232	SN 45233	SN 45234	

5002

Compression Stockings Ccl.2

VENOSAN® 5002
Ccl.II
23–32 mm Hg
(3,1–4,3 kPa)



ARTICLE	COLOUR	LENGTH	TOE	SELF SUPPORTING TOP	S	M	L	XL	
AD	mexico	long	open		SG 46001	SG 46002	SG 46003	SG 46004	
		short			SG 46021	SG 46022	SG 46023	SG 46024	
		long	closed		SG 46011	SG 46012	SG 46013	SG 46014	
		short			SG 46031	SG 46032	SG 46033	SG 46034	
	marokko	long	open		SM 46001	SM 46002	SM 46003	SM 46004	
		short			SM 46021	SM 46022	SM 46023	SM 46024	
		long	closed		SM 46011	SM 46012	SM 46013	SM 46014	
		short			SM 46031	SM 46032	SM 46033	SM 46034	
	sand	long	open		SD 46001	SD 46002	SD 46003	SD 46004	
		short			SD 46021	SD 46022	SD 46023	SD 46024	
		long	closed		SD 46011	SD 46012	SD 46013	SD 46014	
		short			SD 46031	SD 46032	SD 46033	SD 46034	
	black	long	open		SN 46001	SN 46002	SN 46003	SN 46004	
		short			SN 46021	SN 46022	SN 46023	SN 46024	
		long	closed		SN 46011	SN 46012	SN 46013	SN 46014	
		short			SN 46031	SN 46032	SN 46033	SN 46034	
	ADH	mexico	long	open	plain	SG 46041	SG 46042	SG 46043	SG 46044
			short			SG 46061	SG 46062	SG 46063	SG 46064
			long	closed		SG 46051	SG 46052	SG 46053	SG 46054
			short			SG 46071	SG 46072	SG 46073	SG 46074
		marokko	long	open	plain	SM 46041	SM 46042	SM 46043	SM 46044
			short			SM 46061	SM 46062	SM 46063	SM 46064
			long	closed		SM 46051	SM 46052	SM 46053	SM 46054
			short			SM 46071	SM 46072	SM 46073	SM 46074
black		long	open	plain	SN 46041	SN 46042	SN 46043	SN 46044	
		short			SN 46061	SN 46062	SN 46063	SN 46064	
		long	closed		SN 46051	SN 46052	SN 46053	SN 46054	
		short			SN 46071	SN 46072	SN 46073	SN 46074	

ARTICLE	COLOUR	LENGTH	TOE	SELF SUPPORTING TOP	S	M	L	XL
AGH	mexico	long	open	plain	SG 46611	SG 46612	SG 46613	SG 46614
		short			SG 46641	SG 46642	SG 46643	SG 46644
		long	closed	plain	SG 46621	SG 46622	SG 46623	SG 46624
		short			SG 46651	SG 46652	SG 46653	SG 46654
		long		lace	SG 46601	SG 46602	SG 46603	SG 46604
		short			SG 46631	SG 46632	SG 46633	SG 46634
	marokko	long	open	plain	SM 46611	SM 46612	SM 46613	SM 46614
		short			SM 46641	SM 46642	SM 46643	SM 46644
		long	closed	plain	SM 46621	SM 46622	SM 46623	SM 46624
		short			SM 46651	SM 46652	SM 46653	SM 46654
		long		lace	SM 46601	SM 46602	SM 46603	SM 46604
		short			SM 46631	SM 46632	SM 46633	SM 46634
	sand	long	open	plain	SD 46611	SD 46612	SD 46613	SD 46614
		short			SD 46641	SD 46642	SD 46643	SD 46644
		long	closed	plain	SD 46621	SD 46622	SD 46623	SD 46624
		short			SD 46651	SD 46652	SD 46653	SD 46654
		long		lace	SD 46601	SD 46602	SD 46603	SD 46604
		short			SD 46631	SD 46632	SD 46633	SD 46634
	black	long	open	plain	SN 46611	SN 46612	SN 46613	SN 46614
		short			SN 46641	SN 46642	SN 46643	SN 46644
		long	closed	plain	SN 46621	SN 46622	SN 46623	SN 46624
		short			SN 46651	SN 46652	SN 46653	SN 46654
		long		lace	SN 46601	SN 46602	SN 46603	SN 46604
		short			SN 46631	SN 46632	SN 46633	SN 46634
AGG left	mexico	long	open	SG 46151	SG 46152	SG 46153	SG 46154	
AGG right		long		SG 46171	SG 46172	SG 46173	SG 46174	
AT	mexico	long	closed		SG 46221	SG 46222	SG 46223	SG 46224
		short			SG 46231	SG 46232	SG 46233	SG 46234
	marokko	long	closed		SM 46221	SM 46222	SM 46223	SM 46224
		short			SM 46231	SM 46232	SM 46233	SM 46234
	sand	long	closed		SD 46221	SD 46222	SD 46223	SD 46224
		short			SD 46231	SD 46232	SD 46233	SD 46234
	black	long	closed		SN 46221	SN 46222	SN 46223	SN 46224
		short			SN 46231	SN 46232	SN 46233	SN 46234

5003

Compression Stockings Ccl.3

VENOSAN® 5003
Ccl.III
34-46 mm Hg
(4,5-6,1 kPa)



ARTICLE	COLOUR	LENGTH	TOE	S	M	L	XL
AD	mexico	long	open	SG 47001	SG 47002	SG 47003	SG 47004
			closed	SG 47011	SG 47012	SG 47013	SG 47014
	marokko	long	open	SM 47001	SM 47002	SM 47003	SM 47004
			closed	SM 47011	SM 47012	SM 47013	SM 47014
	black	long	open	SN 47001	SN 47002	SN 47003	SN 47004
			closed	SN 47011	SN 47012	SN 47013	SN 47014



VENOSAN®

5000

Compression Stockings

THE MAGIC OF SILVER



VENOSAN
WYDZIAŁ

VENOSAN 5000



VENOSAN®
Compression Stockings

www.venosan.com

VENOSAN® is a registered
trademark of SWISSLASTIC AG
ST. GALLEN, Switzerland

an  Company



SWISSLASTIC AG ST. GALLEN
Rorschacher Strasse 304
CH-9016 St. Gallen
Switzerland