

 chamber

Nebuliser

anywhere
and anytime



- Portable handheld nebuliser
- Lightweight & ideal for home & hospital use
- Quiet in operation
- Detachable parts for easy cleaning
- Built in long life lithium batters with micro USB charging port
- MMAD of 2.5 μ m - 4.5 μ m
- Nebulisation rate of 0.3ml/min - 0.7ml/min

The **e-chamber Portable Nebuliser** is a low cost, lightweight and compact device that is able to provide delivery of asthma and COPD medication, in instances where an asthma spacer or pMDI is not suitable or available. The special technology employed in the **e-chamber nebuliser** allows for smaller medication particles to penetrate the deepest parts of the lungs.



BIRD HEALTHCARE

Exclusive to Bird Healthcare

www.birdhealthcare.com



echamber Nebuliser
anywhere and anytime

Exclusive to Bird Healthcare

Proudly Supporting



www.birdhealthcare.com