



# PULSE OXIMETER

## Accurately measure SpO2 and pulse rate

The Suresense Pulse Oximeter is a noninvasive spot-check measurement of SpO2 and pulse rate. This product is suitable for use in family, hospital (including clinical use in internist/surgery, anesthesia, pediatrics, intensive care etc.), oxygen bar, social medical organisations and physical care in sports.

### Product Details

- Lightweight and Easy-To-Use.
- Colour OLED display, simultaneous display for testing value and plethysmogram
- Adjust the parameters in friendly menu.
- Big fonts mode is convenient for user reading the results.
- Manually adjust the direction of interface.
- Low Battery voltage indicator
- Visual alarm function. Real-time spot-checks
- Automatically switches off when no signal.
- Standard two AAA 1.5V Alkaline battery is available for power supply.
- Advanced DSP algorithm to minimise motion artefacts and improve accuracy of low perfusion SpO2
- Measurement Range: 35-100%
- Resolution: 1%
- Measurement Accuracy: 2% (70%-100%)
- Alarm Range: 50% ~100 %
- Pulse Rate:
- PR Measurement Range: 35-250 BPM
- PR Resolution: 1 BPM
- PR Measurement Accuracy: 2 BPM
- Alarm Range: 30% ~250 BPM
- WARRANTY: 1 Year



Dual colour OLED display



Shock-proof



Low perfusion 0.3%



1300 365 561  
sales@birdhealthcare.com  
www.birdhealthcare.com



BIRD HEALTHCARE